




QuikRead go[®] CRP & CRP+Hb




1 Remove the foil cover from the cuvette. Do not touch the clear flat surfaces in the lower part of the cuvette.



2 Fill the blue striped capillary with the sample (20 µl) up to the white stopper. Ensure that there are no air bubbles in the capillary. Wipe away any residual sample from the outside of the capillary.




3 Place the sample (20 µl) into the buffer solution in the cuvette and dispense by pressing down the plunger. Make sure that the capillary is completely empty.



4 Close the cuvette tightly with the CRP reagent cap. Do not press down the blue inner part of the reagent cap. Perform the test within 2 hours.



5 Choose **Measure** on the display of the QuikRead go instrument.



6 Insert the cuvette in to the measurement well of the instrument. The barcode should be facing you.



7 The result appears on the display and the cuvette rises up automatically when the measurement is completed. The QuikRead go CRP+Hb test shows the Hb result below the CRP result.

Please consult the instructions for use before performing the test.

QuikRead go[®] CRP & CRP+Hb



Kit

Unopened kit at 2–8°C: until the **expiry date** of the kit



Reagent cap tube

Storage at 2–8°C: until the **expiry date** of the kit
Storage at 18–25°C: **1 month** (3 months when 7.5 hours/day)



Cuvettes

Unopened foil pouch at 2–25°C: until the **expiry date** of the kit
Unopened cuvettes after opening the foil pouch:
at 2–8°C: **6 months** / at 18–25°C: **3 months**



Opened cuvette

Let the buffer warm up to 18–25°C before testing.
Perform the test within **2 hours** from opening the cuvette.
The sample is stable in the buffer for **2 hours**.



CRP Controls

Unopened vial at 2–8°C: until the **expiry date**
Opened vial at 2–8°C: up to **2 months**



Hb Control

Unopened vial at 2–8°C: until the **expiry date**
Opened vial at 2–8°C: up to **1 month**

Contact info

Aidian Oy

aidian@aidian.eu

+358 10 309 3000

For orders: orders@aidian.eu

www.aidian.eu

www.quikread.com

Product	Cat. no.
QuikRead go CRP with capillaries and plungers	135171
QuikRead go CRP+Hb with capillaries and plungers	140068
QuikRead go CRP Control	153764
QuikRead go CRP Control High	153763
QuikRead go Hb Control	141154
QuikRead Capillaries 20 µl	67962
QuikRead Plungers	67966
QuikRead go Instrument	133893