



Uricult® Model Chart



Uricult® Products

For diagnosis of urinary tract infections by demonstrating microbes in urine

Uricult 67404

- » CLED medium for total bacterial counts
- » Selective MacConkey medium for gram-negative bacteria



Uricult Plus 67465

- » CLED medium for total bacterial counts
- » Selective MacConkey medium for gram-negative bacteria
- » Selective enterococcus medium for enterococci, group B streptococci may also grow on this medium



Uricult Trio 68197

- » CLED medium for total bacterial counts
- » Selective MacConkey medium for gram-negative bacteria
- » E. coli medium for gram-negative, β -glucuronidase -producing organisms, such as *E. coli*. Certain yeasts may also grow on this medium.



Product list	Cat. No.
Uricult 10 tests	67404
Uricult Plus 10 tests	67465
Uricult Trio 10 tests	68197
Cultura M 230 V	129997 / 05866
Cultura M 110/115 V	132754 / 06266
Cultura M 230 V (with a UK plug)	133192
Cultura M 12 V (with a car plug)	132697
Thermometer, accessory	129981 / 05724
Multirack, accessory	139248

IT'S SO EASY TO USE!



Dip or pour



Incubate



Interpret

For more detailed information, please consult the instructions for use.

Orion Diagnostica develops, manufactures and markets diagnostic test systems for healthcare professionals with over 40 years of experience. Building well-being is the core of our mission.

Uricult® is a registered trademark of Orion Diagnostica Oy.



Orion Diagnostica Oy, Koivu-Mankkaan tie 6 B, P.O. Box 83, FI-02101 Espoo, Finland
Tel. +358 10 4261, Fax +358 10 426 2794, www.oriondiagnostica.com

207-14EN, 10/2017. International version.

Organisms on media



Colony densities *Escherichia coli* CFU/ml

